

Culture Plates 3M™ Petrifilm™



- Ready to use
- 90 % less material volume
- Validated methods
- Time and energy saving
- Increased productivity and reliability
- Low operating costs
- Fast and reliable evaluation

| Description | PK | Cat. No. | Price | Now |
|--|------|-----------|-------|-------|
| | | | EUR | EUR |
| Total bacteria count aerobe bacteria | 100 | 6.237 738 | 99,99 | 88,88 |
| Total bacteria count aerobe bacteria | 1000 | 6.236 670 | 99,99 | 88,88 |
| Rapid test coliforms | 50 | 6.245 052 | 99,99 | 88,88 |
| Total bacteria count coliforms | 50 | 6.269 820 | 99,99 | 88,88 |
| Total bacteria count coliforms | 1000 | 6.235 875 | 99,99 | 88,88 |
| Bacteria count <i>E.coli</i> | 50 | 6.269 819 | 99,99 | 88,88 |
| Bacteria count <i>E.coli</i> and coliforms | 50 | 6.257 512 | 99,99 | 88,88 |
| Bacteria count <i>E.coli</i> and coliforms | 500 | 6.257 513 | 99,99 | 88,88 |
| Bacteria count Enterobacteriaceae | 50 | 6.251 171 | 99,99 | 88,88 |
| Bacteria count Enterobacteriaceae | 1000 | 6.257 509 | 99,99 | 88,88 |
| Bacteria count <i>S. aureus</i> | 20 | 7.658 234 | 99,99 | 88,88 |
| Bacteria count <i>S. aureus</i> | 50 | 6.257 510 | 99,99 | 88,88 |
| Bacteria count <i>S. aureus</i> | 500 | 6.257 511 | 99,99 | 88,88 |
| Rapid test Yeast and Fungi | 50 | 6.282 104 | 99,99 | 88,88 |
| Bacteria count Yeast and Fungi | 100 | 7.657 837 | 99,99 | 88,88 |
| Bacteria count Yeast and Fungi | 1000 | 6.257 514 | 99,99 | 88,88 |

PRECIOS SIN IVA

Microbiological rapid tests 3M™ Dipslides



Easy-to-use indicator for the detection of bacteria, yeast and mold in liquids and on surfaces.

- Easy-to-use and easy to interpret
- Shelf life of 9 months from the date of manufacture

| For | PK | Cat. No. | Price | Now |
|--|----|-----------|------------------|--------------|
| | | | EUR | EUR |
| PCA + TTC + Nutrient | 10 | 6.257 805 | 99,99 | 88,88 |
| PCA + VRBG + CA + TVC + Enterobacteriaceae | 10 | 6.258 926 | 99,99 | 88,88 |
| TSA/Rose Bengal + TVC | 10 | 6.273 017 | 99,99 | 88,88 |
| PCA + OGYE + TVC + Yeast | 10 | 6.287 594 | 99,99 | 88,88 |

Allergenic test

Optionally in a plastic or in 96 Wells/Kit

Lateral flow device in a plastic casset, 25 devices individually packed.
One foil bag with a plate of 96 removable antibody coated wells.

| Description | PK | Cat. No. | Price | Now |
|--|----|-----------|------------------|--------------|
| | | | EUR | EUR |
| L25WAL Walnut Rapid LFD, pack of 25 Tests | 1 | 4.665 443 | 99,99 | 88,88 |
| L25HZL Hazelnut Rapid LFD, pack of 25 Tests | 1 | 4.665 444 | 99,99 | 88,88 |
| L25PST Pistachio Rapid LFD, pack of 25 Tests | 1 | 4.665 445 | 99,99 | 88,88 |
| L25EGG Egg white Rapid LFD, pack of 25 Tests | 1 | 4.665 446 | 99,99 | 88,88 |
| L25GLU Gluten Rapid LFD, pack of 25 Tests | 1 | 4.665 447 | 99,99 | 88,88 |
| L25ALM Almond Rapid LFD, pack of 25 Tests | 1 | 4.665 448 | 99,99 | 88,88 |
| L25MLK Milk Rapid LFD, pack of 25 Tests | 1 | 4.665 449 | 99,99 | 88,88 |
| L25PNT Peanut Rapid LFD, pack of 25 Tests | 1 | 4.665 450 | 99,99 | 88,88 |
| L25SOY Soy Rapid LFD, pack of 25 Tests | 1 | 4.665 451 | 99,99 | 88,88 |
| L25COC Coconut Rapid LFD, pack of 25 Tests | 1 | 4.665 452 | 99,99 | 88,88 |
| L25FSH Fish Rapid LFD, pack of 25 Tests | 1 | 4.665 453 | 99,99 | 88,88 |
| L25CHW Cashew Rapid LFD, pack of 25 Tests | 1 | 4.665 454 | 99,99 | 88,88 |

Sponge-Stick



Sponges and swabs designed specifically for microbial surface sampling. 3M™ Sponge-Sticks are available in wet or dry formats, with and without gloves, to fit many applications. The 1.5" x 3" 3M Sponge-Stick allows sampling without directly handling the sponge.

Sponges and swabs designed specifically for microbial surface sampling.

| Type | Package contents | PK | Cat. No. | Price | Now |
|--|------------------|-----|-----------|------------------|--------------|
| | | | | EUR | EUR |
| Sponge stick, 1.5 x 3 sponge with 10 ml buffered peptone water broth | 4 bags of 25 | 100 | 6.238 511 | 99,99 | 88,88 |
| Sponge stick moistened with 10 ml neutralization buffer | 1 bag of 100 | 100 | 6.266 542 | 99,99 | 88,88 |
| Sponge stick | 5 bags of 20 | 100 | 7.970 470 | 99,99 | 88,88 |

Luminometer Clean-Trace™ LM1



Simple, quick and reliable hand held diagnostic tool to verify cleaning effectiveness and capture, store and manage data for later analysis.

- Compact and easy-to-use
- Rapid results allow immediate corrective action
- To be used with the dry swabs 3M™ Clean-Trace™ UXL100, AQT200 and AQT100

| Type | PK | Cat. No. | Price | Now |
|------|----|-----------|------------------|--------------|
| | | | EUR | EUR |
| | 1 | 4.665 172 | 99,99 | 88,88 |

Dry swabs for Luminometer 3M™ Clean-Trace™ NG3/LM1



6.261 832

| Type | Description | PK | Cat. No. | Price | Now |
|--------------------------------|------------------------|-----|-----------|------------------|--------------|
| | | | | EUR | EUR |
| Clean-Trace™ ATP UXL100 | Dry swabs for surfaces | 100 | 6.243 230 | 99,99 | 88,88 |
| Clean-Trace™ Water AQT100 | Dry swabs, free ATP | 100 | 6.261 832 | 99,99 | 88,88 |
| Clean-Trace™ Water Plus AQT200 | Dry swabs, total ATP | 100 | 6.261 831 | 99,99 | 88,88 |

Soft Carry Case for Luminometer Clean-Trace™ NG3

| Type | PK | Cat. No. | Price | Now |
|--------------------|----|-----------|------------------|--------------|
| | | | EUR | EUR |
| Clean-Trace™ NGSB1 | 1 | 6.286 192 | 99,99 | 88,88 |

Dry swabs Clean-Trace™



6.255 482

| Type | Description | PK | Cat. No. | Price | Now |
|------------|---------------------------------|----|-----------|------------------|--------------|
| | | | | EUR | EUR |
| ALLTECG0 | Surface Protein Test (Allergen) | 60 | 6.261 368 | 99,99 | 88,88 |
| ATP10MLS | Positive control for surface | 10 | 6.251 968 | 99,99 | 88,88 |
| MPRO50 NEU | Surface Protein Plus | 50 | 4.665 741 | 99,99 | 88,88 |

LABORTOPS



Biological Indicators Attest™ for control of steam autoclaves



Biological control involves the use of indicators, which check whether all microorganisms were killed during the sterilisation process. The bio-indicators are unique, as they directly indicate how many microbial spores have actually been killed. Thanks to the reliability and quality of the information, the regular use of a bio-indicator can lead to an improvement in the performance of the sterilisation process and thus to better results.

Features:

- No external laboratory is required for performance evaluation
- Easy handling
- Reasonable to use after repairs or maintenance of autoclaves

Application:

- Program 3.5 minutes, holding time at 121 °C or 134 °C
- After sterilisation, break the vial in order to bring the nutrient solution into contact with the spore strip
- Results in 48 hours
- Meets EN ISO 11138 standard

| Type | PK | Cat. No. | Price | Now |
|--|-----|-----------|----------------------|------------------|
| Biological Indicators Attest™ for steam autoclaves, blue cap | 100 | 6.280 367 | EUR 99,99 | EUR 88,88 |

Multiparameter Indicator 3M™ Comply™



7.970 699

Class 4 indicators with endpoint reaction. With chemical coating, change their colour when the sterilization cycle is successfully completed. Indicator strips for steam sterilization or ethylene oxide sterilization.

- Steam sterilization colour change from yellow to black
- EO sterilization colour change from red to green
- Economical due to the possibility of cutting in half
- Complies with EN ISO 11140-1:2014

| Description | PK | Cat. No. | Price | Now |
|--|-----|-----------|----------------------|------------------|
| Comply™ 1250 chemical indicator strips | 240 | 7.970 699 | EUR 99,99 | EUR 88,88 |
| Comply™ 1251 EO Chemical Indicator Strip | 240 | 6.274 555 | EUR 99,99 | EUR 88,88 |